

**QUANTITATIVE ANALYSIS OF DIFFERENTIAL PHOSPHORYLATION  
IN A *PKC1* OVEREXPRESSION  
STRAIN OF *SACCHAROMYCES CEREVISIAE***

**M<sup>a</sup> Luisa Hernández<sup>1</sup>, Rasmus Hansen<sup>1</sup>, María Jiménez<sup>2</sup>, Humberto Martín<sup>2</sup>,  
Victor Jiménez-Cid<sup>2</sup>, Concha Gil<sup>2</sup> and María Molina<sup>2</sup>**

<sup>1</sup>Unidad de Proteómica, Universidad Complutense de Madrid.  
Parque Científico de Madrid (UCM-PCM).

<sup>2</sup>Departamento de Microbiología II. Universidad Complutense de Madrid

In *S. cerevisiae*, protein kinase C (Pkc1p) is involved in the control of actin polarization and morphogenesis. Pkc1p acts upstream of the cell integrity MAPK pathway. A protein kinase C overexpression strain of *S. cerevisiae* was investigated for differential protein phosphorylation as compared to an isogenic wild type strain.

We have used a phosphoproteomic approach based on quantitative mass spectrometry based on stable isotope labeling with amino acids in cell culture (SILAC).

The *PKC1* overexpression strain was labeled by growth in media containing stable isotopic amino acids, i.e C13 – arginine and C13-lysine, to do differential analysis in a 1:1 protein mixture of both strains using mass spectrometry.

Several phosphopeptide enrichment techniques have been used, and all fractions were analysed by nano – HPLC-MS/MS and neutral loss dependent MS3 on a LTQ mass spectrometer that allowed identification of phosphopeptides using Mascot scoring and quantification with MSquant, a freely distributed program for SILAC quantification.

Of 299 non-redundant phosphopeptides identified and quantified, 93 were upregulated more than 2-fold (average ratio).

The proteomic work was done at the Proteomics Facility of UCM-PCM, a member of ProteoRed Network